

WOUND CARE PRODUCTS

CALCIUM ALGINATE WOUND DRESSING

SILVER CALCIUM ALGINATE DRESSING



SPECIFICATIONS

Algicare Ag utilizes a mix of calcium alginate with the addition of an ionic silver complex that releases silver ions in the presence of wound exudate. As exudate is absorbed, the dressing forms a soft, cohesive gel that intimately conforms to the wound surface, maintaining an optimal moist wound-healing environment.

Advantages:

- High Absorption Capacity
- Broad Spectrum Antimicrobial Effect
- Haemostatic Effect
- Optimum Healing Environment
- Upon contact with wound exudates, sodium alginate rapidly gels which overlays the wound and provides a microenvironment that promotes healing and the formation of granulation tissue.



Specification:

Item No.	Size	Piece/Box
0102521	5×5cm	10/box
0102522	10×10cm	10/box

CALCIUM ALGINATE WOUND DRESSING


SPECIFICATIONS

Algicare TM calcium alginate dressings are derived from brown seaweed. The products are composed of soft, nonwoven fibers shaped as ropes (twisted fibers) or pads/strips (fibrous mats).

Algicare TM calcium alginate dressings are absorbent and conform to the shape of a wound. When packed, an alginate interacts with wound exudate to form a soft gel that maintains a moist healing environment. The release of calcium ion can activate thrombinogen and improve hemostasis.

Advantages:

- H-high absorption capacity of nonwoven fibers and sodium alginate powder.
- Form a gel within the wound to maintain a moist healing environment and over the activity and breed of bacterial
- The release of calcium ion can activate thrombinogen and improve emostasis.
- Facilitate autolytic debridement
- Are easy to apply and remove.



Indications:

- Wounds at risk of infection (by creating an effective barrier to external contaminants).
- Partial-thickness (second-degree) burns
- Diabetic foot ulcers, leg ulcers and pressure ulcer/sores
- Surgical wounds left to heal by secondary intent.
- Traumatic wounds
- Wounds that are prone to minor bleeding, such as those that have been mechanically or surgically debrided.
- Oncology wounds with exudate.

Specification:

Item No.	Size	Piece/Box
0102501	5×5cm	10/box
0102502	10×10cm	10/box
0102503	10×15cm	10/box
0102504	10×20cm	10/box
0102505	15×15cm	10/box
0102506	20×20cm	10/box
0102507	2g×30cm Rope	10/box
0102508	2g×30cm Strip	10/box

CALCIUM ALGINATE WOUND DRESSING

RELIEVE BURN GEL/DRESSING



SPECIFICATIONS

Burncare™ Relieve burn gel/dressing provides superior burn treatment to halt the progress of burn injuries. One of its big advantages over other burn treatments is the dressings don't stick to the burn. Burncare Dressings are sterile and utilize a foam pad that holds the gel and keeps it in place at the wound site, and the open celled foam used in the sterile dressings, to cool the injury effectively. Burncare carries trace amounts of a natural extract called Melaleuca Alternifolia, or Tea Tree Oil. This extract is considered in many countries to have natural anti-germinal features, which help to prevent infection and promote healing.

Advantages:

It also completely non-adhesive, nontoxic, and won't irritate skin. The direct effective is lower the burning heat of burn, scald and sunburn. Less pain. It will provide protection to avoid the damage of burn to healthy organization. Absorb and disperse heat immediately. It is safe to use on first and second degree burns.

Indications:

Use on burns, scalds and sunburn.
Safe to use on first and second degree burns.
Attention: Store at shady and cool place.



Specification:

Item No.	Size	Piece/Box
0102301	3.5g	100/Box
0102302	10×10cm	10/Box
0102303	20×20cm	10/Box
0102304	30×40cm face mask	1/Bag
0102305	55×40cm	1/Bag
0102306	50 ml bottle/spray	10/Box
0102307	120 ml bottle/spray	10/Box
0102308	25g/50g tube	10/Box

PARAFFIN/VASELINE GAUZE DRESSING

SPECIFICATIONS

Paraffin Gauze Dressing / Xeroform Gauze Dressing / Chlorhexidine Gauze Dressing / Composition Biogaze Gauze Dressing

Paraffin Gauze Dressing consists of fabric of Leno weave of cotton, which has been impregnated with white soft paraffin are used to dress denuded areas. It is non-adherent, non allergenic, gamma sterile Dressing. It is soothing, enhances healing in burns, skin grafts skin loss and lacerated wounds.

Advantages:

It is non-adherent, non allergenic. A non medicated tulle dressing effectively supports all phase of wound healing. Impregnated with paraffin. Forms a separate layer between wound and gauze. Facilitates air circulation and hastens healing. Sterilized by Gamma Radiation.

Indications:

For external use only.
Store in a cool place.



Specification:

Item No.	Size	Piece/Box
0102401	5×5cm	100/Box
0102402	10×10cm	10/Box
0102403	10×40cm	10/Box
0102404	20×20cm	10/Box

CALCIUM ALGINATE WOUND DRESSING

HYDROCOLLOID WOUND DRESSING



SPECIFICATIONS

Hydrocolloid wound dressing is made of a wound contact layer of hydrocolloids and a top layer of either semipermeable polyurethane film or polyurethane foam. On contact with wound surface, the hydrocolloid wound dressings create a moist environment, which can reduce pain, promote the formation of granulating tissue, prohibit bacteria infection and speed up the healing process.



Advantages:

Good absorbercy for exudation from the wound.	Easy to apply and remove, avoiding secondary damage to the wound.
Keep the wound moist and accelerate the wound healing, reduce the pain and frequency of dressing change.	Semi-transparent nature makes it easy to observe progress in wound healing.
Waterproof, permeable and prevent the wound from bacteria outside.	Various different sizes for different wounds on different body locations.
Excellent elasticity and self-adhesive, feels comfortable and easy for application.	
The color of dressing can indicate the time for changing a new one.	

Indications:

Can be used to control skin ulcer ,leg ulcer and pressure sores;
Light abrasion wounds;
Second degree burn wounds;
Necrotic wound;
The bordered and standard products are mainly used on light medium exuding pressure sores and leg ulcers;
The thin products are mainly used on dry to light exuding superficial wounds, surgical post operative wounds and abrasions.
It is also used on small wounds towards the end of healing phase.

Specification:

Item No.	Size	Piece/Box
0206101	5×5cm	10/Box
020610	10×10cm	10/Box
0206103	10×15cm	10/Box
0206104	10×20cm	10/Box
0206105	15×15cm	10/Box
0206106	20×20cm	10/Box

HYDROGEL DRESSING

SPECIFICATIONS

Hydrogel dressing is a kind of high polymer material which is similar to life organization .It is net structure with the linking of polymers. It has good biological compatibility. It will not affect the metabolize processing.And the metabolin can discharge by hydrogel. In this moistenvironment,the wound will not have scab.It will heal more quick than in the dry environment.

Applications:

Ulcer	Outside wound	Burns and scald	Scoring	bedsore



Indications:

Not adhere to the wound.
No damage to new organizations.
Kill all kinds of bacterial and avoid to have infection to the wound.

Specification:

Item No.	Size	Piece/Box
0102201	5×5cm	10/Box
0102202	10×10cm	10/Box

CALCIUM ALGINATE WOUND DRESSING

CHITOSAN WOUND DRESSING



SPECIFICATIONS

Chitosan wound dressing is a "green products" which made from the shell of shrimp, crab, etc. The shell made out the disaccharide biological poly-fiber with basicity and high chemical active, so it has good biological consistent, absorbent, ventilate and anti-bacteria etc. The sutures which made from the material, is not necessary to remove the suture. Shell poly-disaccharide fiber dressing can be anti-bacterial, adjust body's function and improve the wound healing, It also use for the burn, the dressing formed a protection lay which control the unnecessary cell The dressing stimulate the skin growing suitable, and minus sear Besides it can solve the problem of the dressing pasting with the wound, it is only necessary to change the outside gauze sponge. So the shell poly-disaccharide medical dress called "artificial skin". Shell poly-disaccharide medical dressing (made by Roosin Medical Co.) proved the wonderful result by clinical test. Some case of operation wound could not be healing used the other dress, but after use Shell poly -disaccharide medical dressing, the wound healed soon. The dressing solved the problem of doctors and the patient pain Now Roosin Medical regards the shell poly-disaccharide medical dressing as one of important projects and produce it in large quantity to solve the patients' pain. Shell poly -disaccharide medical dressing content water 76% tension strength 6.5Mpa, extension degree 210%, tearing strength 5.3KN/M. The dressings are with good biological consistent and anti-bacterial. It can improve the wound healing. It do not paste with wound using this dressing. It is suitable for kinds of wound, inflammation, infection and burn wound.



Advantages:

Sterile package for easy use.
Creating a humid situation, to soft and, to get rid of the dead cell of the wound.
Strong absorbent, it can absorb the 15 times of itself weight.
It can be solved and absorbed, less bleeding, painless.
Anti-bacterial and improve the wound healing.
-It is easy to change the dress as the dressing do not paste with the wound.

SILICONE GEL DRESSING

NON-WOVEN ADHERENT WOUND DRESSING

SPECIFICATIONS

Wound dressing uses non-woven on its back, and acrylate glue above it. It is polyethylene net which has strong absorbent and anti-adhesion. Use it after opening, . It is sterilized in each pack.

Specifications

It can cover on oversewing and unoversewing wound, bruise, cut and scoring .Have fixed drain ,three wipe, detaining needle, puncture needle of vein.

Advantages:

Ventilation well. Avoid dipping to the wound.
Low infective gel, good tolerance.
Unique polyethylene net, change without pain.
-he V specification, which is special designed to the wire guide.
Obey the production standard and high quality of European ethicals.



SILICONE GEL DRESSING

MEDICAL ABSORBENT PAD

SPECIFICATIONS

absorbent material is made of viscose fiber, polyester, etc, and has good absorbency and softness, its surface has special polyethylene grid layer can effectively prevent any pain when replace that wound dressings stickfrom wound skin. As dressing, transfusion text-aid, etc dressings product absorption material. considering the water absorptivity, softness and permeability, Not for human skin allergy, stimulate cell toxicity, and can effectively protective wounds from bacterial infection, have strong absorption performance to achieve national hygiene standards, is a new generation of products series dressings ideal absorbing materials.

Indications:

The absorbent pad is used on bandage ,medical dressing and other dressings.

Advantages:

Has good absorbency and softness.
No pain will caused when exchanging the dressing.
Won't cause allergy,skin irritation or cytotoxicity.

Specification:

Item No.	Size	Piece/Box
010201	5×5cm	10/Box
010202	10×10cm	10/Box



MEDICAL CAD ANTI-ADHESION WOUND DRESSING

SPECIFICATIONS

Medical CAD anti-adhesion wound dressing uses the new high-tech research .The research and manufacture are under the guidance of the Second Military Medical University, Shanghai Second Medical University, and other profession. The raw material ethylic acid is provided by USA HCC. We gained the FDA control certificate from USA which is according to the technical demands of USA medical dictionary and British medical dictionary. We also gained the gold award from China Medical Production Exhibition and Eureka gold award from Brussels World Exhibition. And we gained many special awards from Belgium technology department.

This product is soft, non-toxic, non-smell, high absorbent of exudates ,non-adherent to the wound, and sterilization. This product is demonstrated that its safety, reliability and non-infection on clinic. It reduces time of healing on clinic and less pain to patients .It also reduce the labor intensify of health care workers. It is popular in patients and heath care workers because its lower cost.

Advantages:

Use on surgical wound, acute damaged wound, burn, scald, bruise on surface ,bedsore ,ulcer and other clinic care.



MEDICAL CAD ANTI-ADHESION WOUND DRESSING

Advantages:

National inventive patents. Special to use on first wound care.
 Sterilized disposable package. Don't need to sterilize Full specification. Use on different sizes of wound.
 Drainage and absorbent: the fiber in ethylic acid makes our bodies breath freely .Steam and sweat are easy to drainage .Reduce the possibility of infection and soak of skin.
 Bacteriostat: PH 5.2~5.4. It's disadvantage for bacterial to breed in weak acid.
 Sterilization: It has sterilization function through special fabrication processing.
 Chemotaxis of white cell: Acetate fiber has chemotaxis of white cell which can avoid infection to wound.
 Safety: This product is non-toxic, non-infective ,good biological compatibility, no stimulation to skin.
 Softness and adaptability:: This product is soft and breathable. It can adapt the body curve well. The dressing can avoid to adhere to the wound.
 It will not cause pain when removes or exchanges the wound dressing.
 The comparison of this dressing and traditional dressing

	Medical CAD anti-adhesion wound dressing	Traditional dressing
Resistance to bond	Power to remove=0.2N.good anti-adhesion, less pain to patients .Accelerate healing.	Power to remove>1N.Adrere to the wound .More pain to patients. Prolong the healing time.
Water Absorption	High absorbent to water.Higher than cotton 3~5 times.	Common water absorption. Prolong the healing time.
Ventilation	Good ventilation .Make skin breath freely. Keep the wound in good cure environment.	Common ventilation.
Adaptability	CAD is soft, thin, and comfortable. It can adapt the body curve well. Keep your body move as you can .No damage to the wound.	Gauze or Reticulation dressing which is rough, low flexible, and easy to damage the wound.
Bacteriostat	PH5.2~5.4.It has disadvantage to the breed of bacterial.	No bacteriostat.
Sterilization	It is sterilized by special technology.	No sterilization.
Chemotaxis of white cell	Acetate fiber has chemotaxics of white cell. It improves the anti-infection ability.	Don't have this function.
Biological	Low infection. It can degrade biologically .No damage to environment.	Reticulation dressing has PVC on its surface .It can't degrade, so it pollutes the environment.
Convenience	Personal design, reduce the labor intensify of health care.	ꝑroublesome to use.Quantity of work to health care.Is big.
Advance	International advanced medical idea. Development trend to future.	Can be insteaded.

SILICONE GEL DRESSING

ANTI-ADHESION WOUND DRESSING

SPECIFICATIONS

The umbilicus of newborns is easy for bacteria to enter the body. If the umbilicus cannot be protected well, bacteria will enter it and reproduce quickly, which will make it suppurate and develop fatal ichthemia. According to the clinical observation, as the skin of the newborns is delicate, the common gauze fixed by adhesive plaster is easy to sensitive to the skin and the roughness of the gauze is apt to rub the skin. Besides, the common gauze cannot be fixed tightly, and is easy to slip off in the up-down movements of the belly caused by crying.

Product instrument:

The newly developed dressing is made of CA fiber and cotton elastic bandage. CA fiber is soft, non-toxic, insipid, bacteriostatic, absorbent and not irritative to the skin, and it can decrease bleeding and effusion. Cotton elastic bandage has suitable elasticity, which can keep the dressing on the umbilicus better. Indeed, the newborn umbilicus protection CAD anti-adhesion dressing has the advantages of safety, bacteriostat, air permeability, softness, suitable elasticity, no irritation, and no sensitivity etc.



Advantages:

Safety, bacteriostatic and ventilation.









The dressing is soft ,comfortable to skin ,flexible, no damage to skin and comfortable to babies.

No infection and high water absorption.



SAFE DRESSING TRANSPARENT



Model & Specification	Application	Pcs/Box	Boxes/Ctn	
A101(4.4cm x 4.4cm)	For intravenous injection fixation, especially for infants.	100	40	
Model & Specification	Application	Pcs/Box	Boxes/Ctn	
B101(4.4cm x 4.4cm)	For intravenous injection fixation, especially for infants, with little exudates on the pinhole.	75	40	
Model & Specification	Application	Pcs/Box	Boxes/Ctn	
A101U(4.4cm x 4.4cm)	For intravenous injection fixation, especially for infants, with wing or medicine pot on the IV.	100	40	
Model & Specification	Application	Pcs/Box	Boxes/Ctn	
A103(5.0cmx5.7cm)	For intravenous injection fixation, especially for children.	100	40	
Model & Specification	Application	Pcs/Box	Boxes/Ctn	
A101U(4.4cm x 4.4cm)	For peripheral venous catheter fixing (basic model)	100	40	
Model & Specification	Application	Pcs/Box	Boxes/Ctn	
B105(6.0cmx7.0cm)	For peripheral venous catheter fixing (basic model), with little exudates on the pinhole, also used for baby naval sticker.	75	40	
Model & Specification	Application	Pcs/Box	Boxes/Ctn	
A105U(6.0cmx7.0cm)	For IV fixing with wing or medicine pot.	100	40	
Model & Specification	Application	Pcs/Box	Boxes/Ctn	
A111(8.5cmx10.5cm)	For multi-cavity catheter fixing, pulmonary catheter fixing, hemodialysis catheter fixing.	50	24	

SAFE DRESSING TRANSPARENT




Model & Specification	Application	Pcs/Box	Boxes/Ctn
-----------------------	-------------	---------	-----------

A121(10cmx11.5cm)	For carotid artery injection fixation (basic model).	50	24
-------------------	--	----	----



Model & Specification	Application	Pcs/Box	Boxes/Ctn
-----------------------	-------------	---------	-----------

A123(10cmx12cm)	For multi-cavity catheter fixing, pulmonary catheter fixing, Carotid artery fixing, hemodialysis catheter fixing.	100	24
-----------------	---	-----	----




Model & Specification	Application	Pcs/Box	Boxes/Ctn
-----------------------	-------------	---------	-----------

B123(10cmx12cm)	For multi-cavity catheter fixing, pulmonary catheter fixing, Carotid artery fixing, hemodialysis catheter fixing.	100	40
-----------------	---	-----	----




Model & Specification	Application	Pcs/Box	Boxes/Ctn
-----------------------	-------------	---------	-----------

A405(6.0cmx7.0cm)	For peripheral venous catheter fixing (basic model).	100	40
-------------------	--	-----	----



Model & Specification	Application	Pcs/Box	Boxes/Ctn
-----------------------	-------------	---------	-----------

B405(6.0cmx7.0cm)	For peripheral venous catheter fixing (basic model)	100	40
-------------------	---	-----	----




Model & Specification	Application	Pcs/Box	Boxes/Ctn
-----------------------	-------------	---------	-----------

C133(6.0cmx8.0cm)	For IV fixing with wing or medicine pot.	50	48
-------------------	--	----	----



Model & Specification	Application	Pcs/Box	Boxes/Ctn
-----------------------	-------------	---------	-----------

C134(7.2cmx10cm)	For IV fixing with wing or medicine pot.	50	32
------------------	--	----	----



Model & Specification	Application	Pcs/Box	Boxes/Ctn
-----------------------	-------------	---------	-----------

B205(6.0cmx7.0cm)	For peripheral venous catheter fixing basic model) with little exudates on the pinhole, also used for baby navel sticker.	75	40
-------------------	---	----	----



SAFE DRESSING TRANSPARENT



Model & Specification

Application

Pcs/Box


Boxes/Ctn

A205(6.0cmx7.0cm)

For peripheral venous catheter fixing (basic model).



100

40




SAFE DRESSING TRANSPARENT (WITH PAPER FRAME)

Air permeable, waterproof; applied to dressing and fixing for various wounds, protecting minor cut and scrape wound sterile, and also fixing therapy catheter needle



Air permeable tape used for additional records, convenient to record operator's name, date, and pinhole No... Helps to reinforce fixing and easy for management.

Polyurethane material (PU), designed to be ultra-thin, transparent and air permeable, easy to observe the injection site, wavy coating technology make it better permeability, hypoallergenic, can be fully sticking to



Wavy coating technology, avoid discomfort caused by the tension between the film and the skin when body heat changes and also reduce the possibility of damage to the skin.

Characteristics

Wavy coating technology, transparent and air permeable, easy for the wound to breath, prevent bacteria invasion into the wound. Can quickly remove the over exudates or sweat, which makes it easy to observe the wound. Soft, comfortable, and hypoallergenic, can be applicable to every part of body, strong viscosity, can be kept on the wound nearly for 7 days.



Open the sterile package and take out the transparent dressing, peel off paper packing and stick on the skin. Before sticking, make sure the skin disinfectants on the puncture site are clean and dry, and then remove the paper frame.



Observe the puncture site through the dressing. Make sure that I.V point is at the centre of the dressing, and place the incision of the paper frame right onto I.V catheter.



Stick the dressing to the skin naturally, to not pull or drag it. Press the dressing around the puncture site slightly, then peel off the release paper marked with from the incision, at the same time , press the dressing to smooth it.



4



Smooth the dressing. Write down the patient's name, date, bed no. on record tape, and stick it on one side of the dressing to make the dressing easy to remove. The tape also can be used to enhance the fixing of the tube.

Instruction & Caution

- 1 Clean and disinfect the skin according the doctor's operation instruction before using. After the skin is dry, stick " Safe dressing transparent ".
- 2 Insure the transparent film dressing is big enough, the dressing must be at lease 2.5cm wider than wound or puncture spot so that can be better ticking to the dry and healthy skin.
- 3 New transparent film is required when the I.V. catheter needs to be changed; the USCDC suggests change the transparent film every three days hen it is used for I.V. fixation. The transparent dressing can be used for maximum 7 days when it is used for PICC.
- 4 Changing the transparent film immediately when it is broken or dropped off so as to guarantee the barrier and fixing function of the dressing.
- 5 When changing the dressing, please don't drop off the I.V. catheter.
- 6 The exudation may influence the viscosity of the dressing, please change a new one if there is too much exudation.
- 7 Viscidity will be influenced should there is cleanser, bactericide or antibiotic ointment on the skin.
- 8 Don't drag the dressing when putting it on the skin or it would cause unnecessary hurt to the skin.
- 9 Once there is infection, please tear off the film and do necessary treatment. Make suitable medical care and change the dressing frequently at the ame time or stop using.




Product Description

Transparent dressing transparent, waterproof, applicable for the wrapping and fixing of all kinds of wound, aseptic nursing of postoperative wound, acute injury wound, small incision and laceration wound, and also applicable for the intravenous infusion catheter fixation.

Model No	Specification	Applications	Pcs/Box	Boxes/Ctn
B109	7*10	(Applicable for general surgical wounds, hernia resection, and cecum resection)	100	40
B120	10*10	(Applicable for hernia resection, moles and warts resection endoscopic operation)	100	40
B122	10*15	(Applicable for orthopedic, inter vertebral disk, heart surgery and gynecological surgery)	100	40
B123	10*12	(Applicable for PICC and I II degree bed sore)	100	40
B124	10*20	(Applicable for orthopedic, inter vertebral disk, heart surgery and gynecological surgery)	50	40
B125	10*25	(Applicable for orthopedic, inter vertebral disk, heart surgery and gynecological surgery)	50	40
B126	10*30	(Applicable for orthopedic, inter vertebral disc, heart surgery and thoracic surgery)	100	24
B127	10*35z	(Applicable for orthopedic, post heart operation, wound and joint wound of chest and abdomen)	100	24


SAFE DRESSING - NON WOVEN (NON-WOVEN MEDICAL STERILE)



- 1) Take the Safe Dressing - non woven out of the sterile package and peel off the release paper.
- 2) Make the absorbent pad right on the wound and slightly stick it to the wound, then strip off the whole release paper.
- 3) Press the Safe Dressing - non woven around to smooth it.

Applications

Applicable for the wrapping and fixing of all kinds of wound, aseptic nursing of postoperative wound, acute injury wound, small incision and slight bruise, laceration wound, and also applicable for the intravenous infusion catheter fixation.



Model & Specification	Applications	Pcs/ Box	Boxes/ Ctn
C105(6cm * 7cm)	Applied for incision after.	50	40
C120(10cm * 10cm)	Applied for general surgical wounds, hernia resection, and cecum resection, and cecum resection.	25	20
C122(10cm * 15cm)	Applied for orthopaedic Intervertebral disk, heart surgery and gynecological surgery.	25	20
C124(10cm * 20cm)	Applied for Mastectomy, Intervertebral disk, Cesarean section.	25	20
C125(10cm * 25cm)	Applied for orthopaedic, Intervertebral disk, Heart surgery and Cesarean section.	25	12
C126(10cm * 30cm)	Applied for orthopedic, Intervertebral disk, Heart surgery and thoracic surgery.	25	12
C127(10cm * 35cm)	Applied for orthopedic, post heart operation wound and joint wound of chest and abdomen.	25	12
C132(7cm * 9cm)	Applied for eye surgery, elliptic.	25	20

MEDICAL STICK FIXING (SURGICAL TAPE – MEDICAL TRANSFUSION)

SPECIFICATIONS

Safepore (papertape), SafeTrans (transparent tape) and safe silk (silk tape) are available

Applications

Base material: non woven, polyurethane (PU film), PE punch membrane, silk, plain cloth.

High permeability, hypoallergenic.

Can be used for small cut wound and bruise, dressing fixation, catheter fixing, needle and bandage fixing.



Surgical Tape

Material	Specification	Rolls/Box	Boxes/Ctn
Air permeable non woven	1.25cm * 910cm	24	48
	2.5cm * 910cm	12	48
	5cm * 910cm	6	48
Air permeable paper	1cm * 1000cm	10	108
	1cm * 1500cm	10	108
	2cm * 1000cm	5	108
Microporous coating film	1.25cm * 910cm	24	48
	2.5cm * 910cm	12	48
	5cm * 910cm	6	48



Medical Transfusion

Material	Specification	Rolls/Box	Boxes/Ctn	Boxes/Ctn
Hot rolled non woven	3.5cm * 7cm	1	200	50
	3.5cm * 7cm	5	40	50
	3.5cm * 7cm	10	20	50
	3.5cm * 7cm	25	8	50
PE membrane	4cm * 7cm	1	100	100
	4cm * 7cm	5	20	100
	4cm * 7cm	1	100	100
	4cm * 7cm	5	20	100

SAFE CLOSURE

Specifications

It is a good assistant for sewing up the surgical wound by surgical and cosmetic surgery operation.

Caring the surgical wound after taking out stitches, reducing scar formation.


Bonding and repairing the little wound by hurt in daily life, such as the stab wound, the cut wound and so on.

Steps


Firstly cleaning the wound by physiological saline and keep it dry completely.

Sticking on the wound skin closure vertically from the middle of the wound.

Secondly, applying it from the two edges of wound to the middle for avoiding the wound dehiscence.



0.6cm * 3.8cm



1.2cm * 10cm



Specification	Applications	Pcs/Box	Boxes/Ctn
1.2cm * 10cm	Wound of 5 - 10 cm long	100	10
0.6cm * 10cm	Wound of 5 - 10 cm long	100	10
0.6cm * 3.8cm	Wound of 5 cm long	50	20
0.6cm * 7.5cm	Wound of 5 cm long	50	20
0.3cm * 7.5cm	Wound of 5 cm long	50	20

FUNCTIONAL DRESSING

Characteristics

- 1 Good absorbency for exudation from the wound.
- 2 Keep the wound moist and accelerate the wound healing (reduce the pain and avoid the secondary damage to the wound).
- 3 Waterproof, permeable and prevent the wound from bacteria outside.
- 4 Excellent elasticity and self-adhesive, easy for application.
- 5 Change the dressing when it turns white.
- 6 Easy for peeling off and it will not stick to the wound and keep patients

Application

Low or moderate exuded wound.

Crura venous ulcer.

I-II degree bedsores.

Small superficial burn wound.

Post operation wounds.

Donated skin area.

Superficial wound and beauty operation wound.

Period of granulation and epithelialisation of chronic wound.



ALGINATE WOUND DRESSING

CHARACTERISTICS

- 1 Mixture of natural seaweed fibre with calcium ion, has good histo compatibility.
- 2 Gelization white contact with the exudation of wound and well protect and accelerate the wound healing.
- 3 Rapid absorb the exudation from the wound.
- 4 Soft and comfortable for the skin.
- 5 The release of calcium ion can activate thrombinogen and improve hemostasis.
- 6 No adhesive to wound, no residual on the skin.

Application

Acute and chronic wounds and infected wounds.

Clearing the wound exudation and accelerate the granulation growth and help wound healing.

The alginate padding is used in the lacuna antrum wounds, such as nasal cavity, anus and intestines.



FOAM DRESSING

CHARACTERISTICS


- 1 High absorbency, greatly absorb the exudation from the wound.
- 2 Covered by PU film which can be waterproof and prevent the bacteria infection.
- 3 Create a moist environment and accelerate the wound healing.
- 4 Don't stick on wound when changing the dressing, no. pain.
- 5 Soft, comfortable, also can be used as pad.

Application

High exudation wounds, Ex. Ulcers, bedsores, donated skin areas, burn hurt, etc.



WOUND PLASTER SERIES (PLASTER STICK)



NURSING CONSUMABLES

GLOVES

Characteristics

**COTTON KINT GLOVES
LEAD GLOVES**

**SOAP GLOVES NON WOVEN
INNER GLOVES**

**VINTEX GLOVES
VINYL GLOVES**

**SURGICAL GLOVES
MOLED GLOVES**

1. Latex Examination Gloves Powdered and Powder free

Brand : Vintex Malaysia

Available size: X small, Small , Medium , Large , X large

2. Surgical examination gloves

Available size: X small, Small , Medium , Large , X large

3. Nitrile powder free blue gloves

Brand : Safe glove

Available size: X small, Small , Medium , Large , X large

4. High risk glove

Brand : Safe care

Available size: X small, Small , Medium , Large , X large

5. Vinyl gloves (Powdered and powder free)

6. Polyethylene gloves

7. Cotton Knit gloves

8. Cotton gloves



GAUZE

Characteristics

**GAUZE SPONGE STERILE
GAUZE SPONGE NON-WOVEN**

**GAUZE SPONGE NON-STERILE
X-RAY GAUZE STERILE**

**GAUZE SPONGE WOVEN
X-RAY GAUZE NON-STERILE**

GAUZE SPONGE

10 X10 12PLY V/STERILE 171 100/PKT
10 X10 8PLY V/STERILE 19X151
10 X10 12PLY STERILE 171 100/PKT
10 X 10,BLUE,12PLY 171, STERILE
10 X 10,GREEN,12PLY 171, NON STERILE,100/PKT
10 X15 12PLY 171 V/STERILE 100/PKT
10 X15 12PLY 171 STERILE 100/PKT
5 X 5 12PLY V/STERILE 171 100/PKT
5 X 5 16PLY 171 STERILE DOU. WRAPP. 5/PKT
5 X 5 8PLY V/STERILE 19x15 100/PKT
7.5 X 7.5 12PLY STERILE 19x15 100/PKT
7.5 X 7.5 12PLY 19x15 V/STERILE 100/PKT
30 X 30 16PLY 171STERILE
45 X 45 12PLY 171 STERILE 2/PKT
45 X 45 16PLY 171 V/STERILE 100/PKT
10 X10 6PLY NON/WOVEN V/STERILE
10 X10 6PLY NON/WOVEN STERILE
5 X 5 6PLY V/W STERILE 100/PKT
5 X 5 6PLY V/W V/STERILE 100/PKT
7.5 X 7.5 3PLY V/W, V/STERILE 100/PKT



X-RAY GAUZE

10X10 12PLY 17T STERILE
10X10 12PLY 17T V/STERILE, 100/PKT
10X15 12PLY 17T V/STERILE, 100/PKT
10X15 12PLY 19x15 STERILE, 100/PKT
30X30 12PLY 17T V/STERILE 5/PKT
30X30 12PLY 17T STERILE 5/PKT
30X30 16PLY 17T V/STERILE CHIP/TIE
30X30 16PLY 17T STERILE CHIP/TIE
45X45 12PLY 17T STERILE 5/PKT
45X45 12PLY 17T V/STERILE
45X45 13PLY 17T V/STERILE
45X45 13PLY 17T STERILE 5/PKT
22.5X22.5 13PLY 17T V/STERILE 5/PKT
22.5X22.5 13PLY 17T STERILE 5/PKT

BANDAGES

Characteristics

COHESIVE ELASTIC BANDAGE 1" ,12/PKT
COHESIVE ELASTIC BANDAGE 2",12/PKT
COHESIVE ELASTIC BANDAGE 3"12/PKT
COHESIVE ELASTIC BANDAGE 4" 12/PKT
COHESIVE ELASTIC BANDAGE 6" 12/PKT

COTTON BANDAGE 3" 12/PKT
COTTON BANDAGE 1" 12/PKT
COTTON BANDAGE 2" 12/PKT
COTTON BANDAGE 4" 12/PKT
COTTON BANDAGE 6"-12/PKT
COTTON BANDAGE 8" 12/PK,40S,43X27

CONFIRMING BANDAGE15CM 12/PKT
CONFIRMING BANDAGE 5CM 12/PKT
CONFIRMING BANDAGE 7.5CM 12/PKT
CONFIRMING BANDAGE10CM ,12/PKT

ELASTIC ADHESIVE BANDAGE 2
ELASTIC ADHESIVE BANDAGE 4"
ELASTIC ADHESIVE BANDAGE 3"

LATEX FREE ELASTIC BANDAGE 10cm
LATEX FREE ELASTIC BANDAGE 15cm
LATEX FREE ELASTIC BANDAGE 5CM
LATEX FREE ELASTIC BANDAGE 7.5cm

COTTON BANDAGE RUBBER ELASTIC BANDAGES



CREEP BANDAGES COHESIVE BANDAGES

RUBBER ELASTIC BANDAGE 2"
RUBBER ELASTIC BANDAGE 3"
RUBBER ELASTIC BANDAGE 4"
RUBBER ELASTIC BANDAGE 6"
RUBBER ELASTIC BANDAGE 8"

SPANDEX CREPE BANDAGE 5CM 12/PKT
SPANDEX CREPE BANDAGE 7.5CM 12/PKT
SPANDEX CREPE BANDAGE 10CM,12/PKT
SPANDEX CREPE BANDAGE 15CM 12/PKT
SPANDEX CREPE BANDAGE 20CM 12/PKT

SPANDEX ELASTIC BANDAGE 2 "12/PAK
SPANDEX ELASTIC BANDAGE 3 "12/PKT
SPANDEX ELASTIC BANDAGE 4"12/PKT
SPANDEX ELASTIC BANDAGE-6 " 12/PKT



BAN AID

Characteristics

LONG 19X75MM 100/BOX
LONG 25X75MM100/BOX
ASSORTED 100/BOX
LONG BLUE 19X75MM 100/BOX
LONG ASSORTED 19X75MM 100/BOX
SPOT 100/BOX
X LARGE PUSHBAN, 500/BOX
BLUE, 19X75MM,BLUE, 50/BOX



FACE MASK

Characteristics

EARLOOP MICROTECH/USA-BLUE,50/PKT
EARLOOP MICROTECH/USA-WHITE 50/PKT
PAPER 2PLY 100/PKT
TIE ON - MICROTECH 50/PKT
FACE MASK WITH SHIELD-MICROTECH 50/PKT



SHOE COVER

Characteristics

NON WOVEN
PLASTIC PP MICROTECH 100/PKT
NON WOVEN 100/PKT



ISOLATION GOWN

Characteristics

Isolation gown with cuff and tie blue color
Isolation gown with cuff and tie SMS
Isolation gown SMS, blue, wrap around with button
Isolation gown yellow
Isolation gown PE coated

BED SHEET

Characteristics

Disposable bed sheet blue non sterile size 140x240 cm
Disposable bed sheet white non sterile size 140x240 cm
Disposable bed sheet white sterile size 140x240 cm
Disposable fitted bed sheet
Stretcher covers SMS
Reusable bed sheet



PILLOW COVER

Characteristics

Disposable non-sterile blue color
Disposable non-sterile white color
Disposable pillow cover sterile white color
Reusable heavy duty, blue PVC
Reusable Light duty, white PVC
Reusable cotton pillow cover



MALE URINAL

Characteristics

MALE URINAL WITH LID - GREEN
MALE URINAL WITH LID - WHITE

KIDNEY BASIN



BASIN, EMESIS, GOLD,
700ML (KIDNEY TRAY PLASTIC)
DISPOSABLE BASIN ,EMESIS ,GOLD ,
500ML (KIDNEY TRAY PLASTIC) REUSABLE

DELUXE FEMALE URINAL URINOIR POUR DAME



WASH BASIN



SITZ BATH



EMESIS BASIN



GRADUATED CONTAINER



BEDPAN



ID BAND

Characteristics

ADULT YELLOW 100PCS/BX
ADULT BLUE 100/BX
ADULT GREEN 100/BX
ADULT PINK 100/BX
ADULT RED 100/BX
ADULT WHITE 100/BX
INFANT BLUE 100/BX
PEDIATRIC BLUE 100/BX
PEDIATRIC PINK 100/BX
PEDIATRIC WHITE 100/BX



BODY BAG AND POSTMARTUM KIT




KIT36x90 ADULT STRAIGHT ZIP
KIT 28X46 PAEDIA U ZIPPER
KIT 18X28 BABY-U ZIPPER
36X90 ADULT UZIPPER

CENTRAL LINE TRAYS



BONE MARROW BIOPSY/ASPIRATION TRAYS



The bone marrow trays offer both Illinois-style Sternal-lain needles and our exclusive, enhanced J-style bone marrow biopsy/aspiration needles.


GENERAL PURPOSE INSTRUMENT TRAYS



I.V. START KITS



INCISION AND DRAINAGE TRAYS




SUTURE REMOVAL KITS



SUTURING TRAY



J CLOTH



31 CMX37CMX70GMS WHITE 500PCS/CTN
 31X37 CM 45GSM BLUE 500/CTN -AL AIN
 30X60CM- 75GSM-BLUE 500/CTN
 31X37CM- 75GSM - WHITE
 30X60CM- 75GSM-GREN 500/CTN
 30X60CM- 75GSM-RED 500/CTN
 30X60CM- 75GSM-WHITE 500/CTN

WASH CLOTH



MOUTH CARE KIT



DISPOSABLES PANTIES

LARGE,6/PKT
 MEDIUM,6/PKT
 NON WOVEN
 X LARGE, 6/PKT
 XXLARGE



CRASH CART LOCK




CRASH CART LOCK BLUE
 CRASH CART LOCK RED WITH NO.
 CRASH CART LOCK WHITE
 CRASH CART LOCK-BLACK
 CRASH CART LOCK YELLOW
 SECURE PULL SEAL YELLOW 14.25"

NURSING CONSUMABLES

BETADINESWAB STICK



ALCOHOL SWAB



LEMON GLYCERIN SWAB STICK



COTTON ROLL& COTTON BALL

COTTON 1KG ROLL
COTTON 500G ROLL
COTTON 250GM
COTTON 50G ROLL
COTTON 100GM ROLL
COTTON BALL 100/PKT
COTTON BALL, STERILE 100/PKT
COTTON BALL, STERILE 20/PKT
COTTON BALL, STERILE 5/PKT
COTTON BALL 200/PKT

DISPOSABLE DOCTOR COAT& LAB COAT

LAB COAT LARGE DISPOSABLE
LAB COAT XLARGE DISPOSABLE
LAB COAT MEDIUM DISPOSABLE
LAB COAT SMALL DISPOSABLE
LAB COAT XXL DISPOSABLE

STRETCHER COVER



STRETCHER FITTED COVER SHEET 100X240

DRAW SHEET



DRAW SHEET LIGHT DUTY, 60" X 42"
DRAW SHEET HEAVY DUTY, 60" X 42"
DRAW SHEET 80 X 175CM 100/CTN

DRESSING PAD



AMNICATOR



AMNICATOR SWAB STIC MEDICAL WARE 100/PK

PROTOSCOPE




PROCTOSCOPE+PUSHER-EOS
PROTOSCOPE LARGE EOS-BLUE
PROTOSCOPE MEDIUM EOS-GREEN
PROTOSCOPE EOS, ORANGE
PROTOSCOPE SMALL EOS-YELLOW
PROTOSCOPE -WITH LIGHT

EYE PAD



EYE PAD 55 MM X 75 MM BIN,50/PKT
EYE PAD 60MM X 80 MM LARGE ,50/PKT
EYE PAD 20/PKT, SELF ADHESIVE

RYLES TUBE



FEEDING TUBE



FEEDING TUBE 10FR-XRAY BLISTER
FEEDING TUBE 12FR-XRAY PE BAG
FEEDING TUBE 14FR-XRAY PE BAG
FEEDING TUBE 16FR-XRAY PE BAG
FEEDING TUBE 18FR-XRAY PE BAG
FEEDING TUBE 20FR-XRAY PE BAG
FEEDING TUBE 22FR-XRAY PE BAG
FEEDING TUBE 6FR-XRAY BLISTER
FEEDING TUBE 7FR
FEEDING TUBE 8FR-XRAY BLISTER

SUCTION CONNECTION TIP



SUCTION CONNECTION TUBES




SUCTION CONNECTION TUBE AA 1.8M
SUCTION CONNECTION TUBE AA 3.6M
SUCTION CONNECTION TUBE AB 3.6M

SUCTION CONNECTION TIP



MEDICINE ZIPPER BAGS SIZE:100 X 120MM
MEDICINE CUP 30ML / MEDICINE CUP 60ML
MEDICINE CUP 30ML INJECTION MOULD
(vacuum foaming)

PEDIATRIC URINE COLLECTOR




URINE COLLECT BAG PEDIATRIC 100ML

NEEDLES



BLUNT NEEDLE 19G X 1.5" / BLUNT NEEDLE 18G X 1.5"
BLUNT NEEDLE WITH FILTER 18G X 1.5"
NEEDLE FOR EAR IRRIGATION, 15G, STERILE, BLUNT.
NEEDLE EXTR SHOR NEEDLE 30X3/8
NEEDLE FREE VALVE SMARTSITE#2000E
NEEDLE HOLDER WORLD DENTO / NEEDLE ZAP

SCALP VEIN SET



SCALP VEIN SET-18G 100/PKT
SCALP VEIN SET-19G 100/PKT
SCALP VEIN SET-20G,100/PKT
SCALP VEIN SET-21G 100/PKT
SCALP VEIN SET-22G,100/PKT
SCALP VEIN SET-23G,100/PKT
SCALP VEIN SET-24G,100/PKT
SCALP VEIN SET-25G,100/PKT
SCALP VEIN SET-26G ,100/PKT
SCALP VEIN SET-27G,100/PKT

SYRINGES



SYRINGE PUMP-SYRINGE PUMP-Perfusor Compact
 SYRINGE 10ML 18G X 1.5 LL,100/BX
 SYRINGE 10ML 18G X 1.5 LS,100/BX)
 SYRINGE.10ML.21GX11.5-LS-2PART,100/BX
 SYRINGE 10ML 21G X 1.5-LL 100/BX
 SYRINGE 10ML 21G X 1.5 -LL 2PART
 SYRINGE 10ML 21G X 1.5 LS 100/BX
 SYRINGE 10ML TERUMO 100/ BOX
 SYRINGE 10ML W/O NEEDLE LL
 SYRINGE 10ML W/O NEEDLE LS,100/BX
 SYRINGE 10MLW/O NEEDLE-LS-AMBER(55PCS/BX)
 SYRINGE 1ML 26G 1/2 LS
 SYRINGE 1ML 27GX0.50 LS,100/BX
 SYRINGE 1ML DISPOVAN ,100/BX
 SYRINGE 1ML 27GX1/2 LS TUBERCULIN100/BX
 SYRINGE 1ML W/O NEEDLE-AMBER.LS,100/BX
 SYRINGE 20ML TERUMO
 SYRINGE 20ML 18GX1.5 LL, 50/BX
 SYRINGE 20ML 18GX1.5 LS,50/BX
 SYRINGE 20ML 21GX1.5 LL 50/BX
 SYRINGE 20ML 21GX1.5 LS,50/BX
 SYRINGE 20ML W/O NEEDLE LL
 SYRINGE 20ML W/O NEEDLE- AMBER W/CAP LS-50/BX
 SYRINGE 20ML W/O NEEDLE CONCENTRIC- LS,50/BX
 SYRINGE 2.5ML 23GX1.25
 SYRINGE 2ML 23GX1.25-2PART.LS,100/BX)
 SYRINGE 2ML CONCENTRIC W/O NEEDLE- LL,100/BX
 SYRINGE 1ML 30GX1/2 INSULIN- MT(100PCS/BX)
 SYRINGE 30ML 18GX1.5 LL
 SYRINGE 30ML 18GX1.5 LS
 SYRINGE 30ML 21GX1.5 ECCENTRIC
 SYRINGE 30ML 21GX1.5 LL,40/BX
 SYRINGE 30ML 21GX1.5 LS,40/BX
 SYRINGE 30ML W/O NEEDLE, LL,30/BOX
 SYRINGE 30ML W/O NEEDLE LS
 SYRINGE 3ML AMBER W/ONEEDLE
 SYRINGE 3ML TERUMO
 SYRINGE 3ML 23GX1.25 LL100/BX
 SYRINGE 3ML LL W/O NEEDLE 100/BX
 SYRINGE 3ML ARTERIAL BLOOD GAS SAMPLER-3303
 SYRINGE 100ML CATHETER TIP
 SYRINGE 3ML LS W/O NEEDLE,100/BX
 SYRINGE 3ML LS W/ONEEDLE .AMBER(110PCS/BX)

SYRINGE 50/60 ML 18GX1.5 LS (50PCS/BX)
 SYRINGE 50/60 ML 21GX1.5 LS,[50/BX]
 SYRINGE 50ML 18GX1.50 LL 25/BOX
 SYRINGE 50ML 18GX1.5 LS
 SYRINGE 50ML 21GX1.5 ECCENTRIC
 SYRINGE 50ML CATHETER TIP 25/BX
 SYRINGE 50ML W/O NEEDLE LS
 SYRINGE 50ML W/O NEEDLE TERUMO
 SYRINGE 50ML W/O NEEDLE LL-50/BX
 SYRINGE 5ML TERUMO
 SYRINGE 5ML 23GX1.25 LL 100/BX
 SYRINGE 5ML 22GX1.5 LL 100/BX
 SYRINGE 5ML 22GX1.5 LS
 SYRINGE 5ML 22X1.5 -2PART.LS(100PCS/BX)
 SYRINGE 5ML W/O NEEDLE LL
 SYRINGE 5ML-AMBER
 SYRINGE 5ML W/O NEEDLE LS CONCENTRIC
 SYRINGE 5ML W/O NEEDLE-LS-AMBER(75PCS/BX)
 SYRINGE 60ML W/O NEEDLE
 SYRINGE -BULB SYRINGE 60 ML
 SYRINGE MT 50/60CC 21GX1.50"LL
 SYRINGE MT50/60CC-CATHETERTIP,25/BX

IV CANNULA




IV CANNULA 14G
 IV CANNULA 16G
 IV CANNULA 17G
 IV CANNULA 18G
 IV CANNULA 18G CATHY
 IV CANNULA BD 20G
 IV CANNULA 20G
 IV CANNULA 22G
 IV CANNULA 22G CATHY
 IV CANNULA 22G SAFETY MEDIFLON
 IV CANNULA 24G
 IV CANNULA 26G
 IV CANNULA 24G CATHY
 IV CANNULA 20G CATHY

BLOOD INFUSION SET



IV SET



IV SET 20 DROPS Y TYPE MCT
IV SET 20 DROPS MCT
IV SET 60 DROPS Y TYPE MCT
IV SET 60 DROPS MICROTECH

IV BURET SET



IV BURETTE SET 100ML M.TECH
IV BURETTE SET 150ML M.TECH

UMBILICAL CORD CLAMP



3 WAY STOP CLAMP



EXTENSION TUBE



TABLET BOTTLE



TABLET BOTTLE AMBER 8DR [30ML] ,410/CTN
TABLET BOTTLE AMBER 13DR (48ML) ,320/CTN
TABLET BOTTLE AMBER 16DR
TABLET BOTTLE AMBER 20DR ,360/CTN
TABLET BOTTLE AMBER 30DR [112ML],240/CTN
TABLET BOTTLE AMBER 40DR 180/CTN
TABLET BOTTLE AMBER 60DR,110/CTN

NURSES CAP



SURGEON CAP - MICROTECH 100/PKT
SURGEON CAP-MEDPLUS 100/PKT
NURSES CAP BLUE 100/PKT,20PKT/CTN
NURSES CAP WHITE 100/PKT-BINALI, 20/PKT

EYE SHIELD



EYE SHEELD - CLEAR UNIVERSAL 50/PKT

PLASTIC APRON



PLASTIC APRON .750X1200MM EP:100/PKT
PLASTIC APRON 800X1300MM EP, 100/PKT
PLASTIC APRON - ASSORTED COLOURS
PLASTIC APRON ROLL BLUE&WHITE

BULB SYSRINGE



EAR BULB SYRINGE 30MLI NON STERILE
EAR BULB SYRINGE 60ML WITH NOZZLE RUBBER NON-STERIL

ARM BOARD WITH STRAP



ARM BOARD BENDABLE MEDIUM 3X6 #1008-50
ARM BOARD LARGE 3X9 #1009-50
ARM BOARD WITH STRAP - ADULT
ARM BOARD BENDABLE SMALL 2X4 #1004-50

SILVER NITRATE STICK



IV SET 20 DROPS Y TYPE MCT
IV SET 20 DROPS MCT
IV SET 60 DROPS Y TYPE MCT
IV SET 60 DROPS MICROTECH

DISPOSIBLE SLIPPER




TONGUE DEPRESSOR



TONGUE DEPRESSOR NON-STERILE 100/PKT
TONGUE DEPRESSOR PEDIA
TONGUE DEPRESSOR PLASTIC , 50/PKT
TONGUE DEPRESSOR STERILE 6"100/PKT

CATHETER MOUNT



ULTRASOUND PROBE COVERS



ULTRASOUND PROBE COVER
ULTRASOUND PROBE COVER STERILE
ULTRASOUND PROBE COVER PC04 (290MM LENGTH), 100PCS/BOX

EAR THERMOMETER PROBE COVERS




ORAL THERMOMETER PROBE COVERS



OTOSCOPE COVER



LAMINATED COUCH ROLL



COUCH ROLL 60 X 50 PLASTIC LAMINATED 6/CTN
COUCH ROLL 50 X 50 PLASTIC LAMINATED 6/CTN
COUCH ROLL 80 X 50 PLASTIC LAMINATED 6/CTN

COUCH ROLL



COUCH ROLL 50 X 50 2PLY 12/CTN AJM
COUCH ROLL 2PLY 12/CTN-UAE
COUCH ROLL 40 X 50 2PLY 9ROLL/CTN
COUCH ROLL 50 X 50 CLSIC 9/CTN
COUCH ROLL 50 X 50 12/CTN
COUCH ROLL 50 X 50 ERREBI ITALY 9/CTN

MAXI ROLL



MAXI ROLL 12/CTN, AL AIN
MAXI ROLL 2PLY 6/CTN
MAXI ROLL DESPENSOR

TOILET ROLL



TOILET ROLL -SILVERCITY
TOILET ROLL-MINITORK, 12/CTN

MAXI ROLL DISPENSERS



C FOLD DISPENSERS



SOAP DISPENSERS



FACE TISSUE



FOAM CUP



PLASTIC CUP



THERMAL PRINTER PAPER



THERMAL PAPER FOR THE ANALYZER FUKUDA SANGYO ST75, 57 X20MM
THERMAL RECORDING PAPER (HD PAPER)-MITSUBISHI

ULTRASOUND PRINTER PAPER



DIAPER




DIAPER X-LARGE WITH INDICATOR
DIAPER LARGE WITH INDICATOR, 80/CTN
DIAPER SMALL WITH INDICATOR, 80/CTN
DIAPER MEDIUM WITH INDICATOR, 80/CTN

UNDER PAD



UNDER PAD 60 X 45 CM 38G
UNDER PAD 60 X 60 CM 60G-240/CTN
UNDER PAD 60 X 90 CM 70G-150/CTN
ER PAD 60 X 90 CM 80G-150/CTN

ECG PRINTER PAPER



ECGPAPERROLL 5 CM # F 1 576
ECGPAPERROLL 60MM-FS722
ECG PAPER- B2-180
ECG PAPER M 2483A (A4 SHEET) 200PCS/PKT
ECGPAPERROLL50X30 X5CM40457C/DBORD06-4 -#13312
ECGPAPERROLL50X30 X6CM 40457C/D BORD06-4 -#21697
ECGPAPER107,75X23 LIFA PACK11RS100 -21885
ECGPAPER210X150MM 200 SHEETS2 FOLD
ECGPAPER 90MM X 70MM X 400 EP-80, SH. CODE: 21773
ECGPAPERROLL210X30CARDIOLINEDELTAEX -13328
ECGPAPER63 x30x 75 FUKUDA OP122TEX#21599
ECGPAPER63 x30x 75 FUKUDA OP122TEX-13145
ECGPAPER N3-735
ECGPAPER N4-879
ECGPAPER210MMX280MMX56m,200S
ECGPAPER210MMX300MMX200SM1709A HP-21710
ECGPAPER80 X 70MM X314 (P8000)#21543
ECGPAPER210x30m BIGSIZEFORHPC14-841

EMERGENCY BLANKET



EMERGENCY BLANKET GOLD-SILVER SURFACE 1600X2100MM
EMERGENCY BLANKET SILVER SURFACE 1422X2100MM

AUTOCLAVE TAPE



AUTOCLAVE TAPE 1"
AUTOCLAVE TAPE 3/4"

CARE BAG



COLD PACK



COLD AND HOT PACK
COLD PACK INSTANT 5"X6" 84/CTN
COLD PACK INSTANT 6" X 9" 48/CTN

URINE BAG



BEDPAN LINER



VOMIT BAG

